INSTRUCTION FOR USE

DEVICE DESCRIPTION

The Intraprostatic Temperature Probe (IP) monitor the temperature inside the prostate during treatment, see CoreTherm® IFU L-2013-009/ CoreTherm® Eagle IFU L-011.

- The Intraprostatic Temperature Probe is a plastic tube containing temperature sensors. The tube has a Luer fitting with an integrated locking nut that secures the probe to the CoreTherm® Catheter.
- The Intraprostatic Temperature Probe contains four temperature sensors.
- Three temperature sensors are located at the tip of the probe, these monitor the temperatures inside the prostate during treatment. A fourth sensor, located further down from the tip monitors the temperature of the water circulating in the catheter during treatment.
- The Intraprostatic Temperature Probe has a colour coded (green) control unit connector. The connector has a hardware key that prevents incorrect attachment of the probe to the control unit. An electronic chip located in the connector contains a unique serial number that identifies the specific temperature probe in use.

The Intraprostatic Temperature Probe is available in two different versions, Short and Standard, see Fig 1. The specific prostate dimensions of the patient determine which version should be used.

Intraprostatic Temperature Probe Short (yellow)

- Prostate lengths between 25-34 mm and/or prostate weights between 20-30 g.

Intraprostatic Temperature Probe Standard (blue)

- Prostate lengths ≥ 35 mm and prostate weight from 50 g.

Intraprostatic Temperature Probe Short and Standard

- Both Short and Standard IP can be used with prostate lengths ≥ 35 mm and prostate weights between 31-49 g.

		Prostate length	
		≥ 35 mm	25-34 mm
Prostate weight	≥ 50 g	Standard (blue)	Short (yellow)
	31-49 g	Standard (blue) Short (yellow)	Short (yellow)
	20-30 g	Short (yellow)	Short (yellow)

Valid prostate dimensions for treatment are prostate weights ≥ 20 g and prostate lengths ≥ 25 mm. Patients with prostate dimensions below these values should not be treated with CoreTherm®.

INTENDED USE

The Intraprostatic Temperature Probe monitors the temperature inside the prostate during treatment with CoreTherm®.


Figure 1. IP Standard (blue) and IP Short (yellow)

CoreTherm® with accessories is intended for treatment of lower urinary tract symptoms (LUTS) due to benign prostatic enlargement (BPE) using microwave thermotherapy by, or under the supervision of, a trained medical practitioner.

Intended patient population

Adult male patients with BPE.

Intended user profile

To be used by, or under supervision of, a trained medical practitioner.

Intended use environment

Hospitals and health care centers.

INDICATIONS

- The Intraprostatic Temperature Probe is used for monitoring the temperature inside the prostate and the temperature of the water circulating in the catheter during CoreTherm® treatment.
- See CoreTherm® IFU, L-2013-009/ CoreTherm® Eagle IFU L-011 for indications for CoreTherm® treatment.

CONTRAINDICATIONS

- No known contraindications with the Intraprostatic Temperature Probe.
- See CoreTherm® IFU L-2013-009/ CoreTherm® Eagle IFU L-011 for contraindications for CoreTherm® treatment

WARNINGS

There is a risk of complications if the Intraprostatic Temperature Probe is not handled carefully and correctly. The risk of complications can be minimised by using the product as described in this instruction for use.

- National requirements for reprocessing must be complied with if stricter than described reprocessing procedure in this IFU.
- The probe should always be visually checked after each step of cleaning/disinfection and sterilisation for cracks, discoloration, unnatural bending, or loose parts.
- A damaged probe shall be disposed.
- The Intraprostatic Temperature Probe must be cleaned, thermally disinfected and sterilised before use.
- It is important to ensure that no residues of cleaning agents remain on the probe before it is inserted into the patient's urethra.
- Insufficient retraction of the Intraprostatic Temperature Probe before the catheter is removed may cause patient injuries. It is therefore important to retract the Intraprostatic Temperature Probe at least 50 mm before the CoreTherm® Catheter is removed.
- Ensure that the Intraprostatic Temperature Probe is fully inserted and locked to the CoreTherm® Catheter to be correctly positioned in the patient's prostate before treatment is started. Inaccurate temperature feedback may result in the application of insufficient or excessive heat to the prostate.
- Not calibrated or incorrectly calibrated temperature probes can result in incorrect temperature measurements that may jeopardise the safety of the patient and cause permanent injury.
- CoreTherm®/SE: Manufacturer's probe calibration instructions given are strictly followed (see CoreTherm IFU L-2013-009).
- CoreTherm® Eagle: Probe calibration shall be done by ProstaLund. Contact ProstaLund for further information.

PRECAUTIONS

- Do not bend the Intraprostatic Temperature Probe to a radius less than 50 mm.
- During cleaning, special attention should be paid to hard-to-reach areas, such as the coupling.
- After cleaning/ disinfection, make sure that no blood residue is left in the blue or yellow connector or in any other areas.

DIRECTIONS FOR USE

Prior to treatment

- Inspect the Intraprostatic Temperature Probe for cracks or any sign of physical damage before use. Exchange any damaged Intraprostatic Temperature Probe for a fully operational, undamaged probe before starting the treatment.
- The Intraprostatic Temperature Probe can be reused and is supplied in non-sterile packaging. Prior to every use the Intraprostatic Temperature Probe shall be cleaned, disinfected and sterilised by autoclave according to the instructions below.

Treatment

- Insert the sterile Intraprostatic Temperature Probe fully into the CoreTherm® Catheter, then withdraw the probe so that its tip is visible but not protruding from the catheter. Do not secure the probe to the catheter yet.


Figure 2. Withdraw the probe

- When the CoreTherm® Catheter is fully inserted, and the balloon is inflated, gently withdraw the CoreTherm® Catheter until the balloon meets the bladder neck and thus is secured in place.
- Hold the CoreTherm® Catheter to straighten it and insert the Intraprostatic Temperature Probe into the prostate. Advise the patient that he may experience a slight pricking sensation.
- After the Intraprostatic Temperature Probe is inserted into the prostate, secure it to the catheter by turning the locking nut of the probe connector in the clockwise direction.


Figure 3. Secure the Intraprostatic Temperature Probe

- Connect the Intraprostatic Temperature Probe to the green socket of the CoreTherm® System, a click sound confirms the connection. For further instructions please read the CoreTherm® IFU L-2013-009 or CoreTherm Eagle IFU L-011.

After treatment

- Release the Intraprostatic Temperature Probe from the CoreTherm® Catheter by turning its locking nut in the counter-clockwise direction. Then carefully withdraw the probe until it entirely resides within the catheter (minimum 50 mm).

- Visually check that the Intraprostatic Temperature Probe is intact after treatment.
- Disconnect the Intraprostatic Temperature Probe from the green socket of the control unit by pulling it outwards.
- The Intraprostatic Temperature Probe shall be cleaned, disinfected, and sterilised prior to every use.

NOTE: The cleaning should be conducted as soon as possible after use of the Intraprostatic Temperature Probe to avoid blood and particles to dry in on the probe.

- The Intraprostatic Temperature Probe can be reused, provided that the instructions given in this section are followed. Since the temperature probe is fragile, it is very important that it is handled with care.

NOTE: Do not bend the Intraprostatic Temperature Probe to a radius less than 50 mm.

AUTOMATED CLEANING AND THERMAL DISINFECTION

1. Remove blood stains and gross soil from the probe using a soft brush with synthetic fibers under running lukewarm tap water. Water quality to be used shall comply with AAIM TIR34. A small pipe brush can be used for the narrow areas.
2. Use only validated washer-disinfector machines with approved efficacy (e. g. CE mark or FDA clearance and validation according to ISO 15883).
3. Use only cleaning agent intended for use in washer-disinfector (e. g. Neodisher® MediClean Forte). Do not exceed the concentration, temperature, and contact time recommended by the detergent manufacturer.
4. Load the probe into the washer-disinfector. Avoid contact between the probe and the other devices as movement during washing could cause damage and washing action could be obstructed.
5. Make sure not to bend the Intraprostatic Temperature Probe to a radius less than 50 mm.
6. Run the washer-disinfector cycle. Water quality to be used shall comply with AAIM TIR34. A thermal disinfection program of 5 minutes in 90°C is required.

7. Upon completion, visually inspect the probe for any damage or remaining soil or blood stains. If soil or blood stain remains, repeat the cleaning process.
8. If additional drying is required, arrange instruments in a clean area or use a soft, lint-free cloth.

NOTE: Chemical disinfection programs are not recommended due to the potential for chemical residues to remain on the probe. These residues could interfere with sterilisation efficacy.

STERILISATION

The Intraprostatic Temperature Probe can be reused and is supplied in non-sterile packaging. The Intraprostatic Temperature Probe shall be cleaned, disinfected and sterilised by autoclave prior every use. The probe shall be sterilised by autoclave according to the instructions below.

Autoclave sterilisation

Dynamic-Air-Removal Steam Sterilisation process follows the ISO 17665-1 standard for steam sterilisation according to following parameters:

1. Place the cleaned probe in a Tyvek® bag.
2. Place the bag with the probe in the Autoclave.
3. Sterilise according to the following cycle:

Temperature:	121°C
Pre-vacuum	
Time:	min 15 minutes (max 30 minutes)
Minimum dimension of packaging:	205x270 mm

CARE AND STORAGE

- The probe must be stored in a clean and dry place, where the temperature does not fall below 5°C or rise above 40°C, and humidity is between 30% and 75% RH.
- Do not bend the Intraprostatic Temperature Probe to a radius of less than 50 mm.

REUSE

The Intraprostatic Temperature Probe is a reusable device that is delivered unsterile. Recommended lifetime use is 10 times. It is not possible to reuse the probe more than 15 times. The Intraprostatic Temperature Probe must be cleaned, disinfected and sterilised after every use.

SAFE DISPOSAL

The Intraprostatic Temperature Probe shall be disposed of according to the local regulations applying to hazardous biological substances.

RESPONSIBILITIES

- Sterilisation process and equipment characterisation (description, specification, ancillary specification, and service definition) should be handled by the Steriliser manufacturer, in collaboration with the developer of the sterilisation process.
- Validation (Installation qualification, operational qualification, performance qualification as well as review and approval of validation) should be handled by the organisation with responsibility of sterilising the product.
- Routine monitoring and control (load, configuration, process monitoring, record generation, periodic testing, and record retention) should be handled by the reprocessing facility.
- Product release from Sterilisation (record review, indicator testing, product disposition and corrective action) should be handled by the reprocessing facility.

Maintaining sterilisation process effectiveness should be handled by the reprocessing facility together with the organisation sterilising the product.

REPORTABILITY

Any serious incident that has occurred in relation to the medical device shall be reported to ProstaLund AB and to the authority having jurisdiction in your locale.

SYMBOLS

Manufacturer

Manufacturing date, YYYY-MM

Type BF applied part

Read instructions for use

Range of humidity

Temperature limits

Article number

Serial number

Medical Device - Indicates the item is a medical device

ORDERING INFORMATION

Device	Article Number
Intraprostatic Temperature Probe, AC Standard (IP Standard autoclavable)	IP807015
Intraprostatic Temperature Probe, AC Short (IP Short autoclavable)	IP807115

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