
INSTRUCTION FOR USE

DEVICE DESCRIPTION

The CoreTherm® Catheter (CA) is a single use device that has separate channels for the Microwave Antenna, the Intraprostatic Temperature Probe, the water inlet and outlet and for inflating the balloon.

The CoreTherm® Catheter is to be used with CoreTherm®, see CoreTherm® IFU L-2023-009 / CoreTherm® Eagle IFU L-011.

- The CoreTherm® Catheter is fitted with a balloon close to its tip.
- The balloon resides inside the patient's bladder and is inflated to keep the catheter anchored at the neck of the bladder during treatment.

The purpose of the CoreTherm® Catheter is to house the Microwave Antenna (AN) and the Intraprostatic Temperature Probe (IP), and to remove heat loss from the Microwave Antenna during treatment.

The catheter, the water tubes and the water reservoir are preassembled and must not be separated.


Figure 1. The CoreTherm Catheter with the water tubes and the water reservoir.

The CoreTherm® Catheter is supplied sterile in a sealed and sterile package with a unique catheter serial number (SN) and is intended for single use.

The CoreTherm® Catheter is available in two different versions, Short and Standard. The two different versions have colour-coded inlet ports for the Intraprostatic Temperature Probe and the Microwave Antenna, see Fig 2. Prostate dimensions measured with transrectal ultrasound determine which version should be used.


Figure 2. Colour-coded inlet ports: Short (yellow) and Standard (blue).

CoreTherm Catheter Short (yellow)

- Can be used with prostate lengths between 25-34 mm and/or prostate weights between 20-30 g.

CoreTherm Catheter Standard (blue)

- Can be used with prostate lengths ≥ 35 mm and prostate weights ≥ 50 g.

CoreTherm Catheter Short and Standard

- Both Short and Standard can be used with prostate lengths ≥ 35 mm and prostate weights between 31-49 g.

Valid prostate dimensions for treatment are prostate weights ≥ 20 g and prostate lengths ≥ 25 mm. Patients with prostate dimensions below these values should not be treated with CoreTherm®.

		Prostate length	
		≥ 35 mm	25-34 mm
Prostate weight	≥ 50 g	Standard (blue)	Short (yellow)
	31-49 g	Standard (blue) Short (yellow)	Short (yellow)
	20-30 g	Short (yellow)	Short (yellow)

INTENDED PURPOSE

Intended use

The CoreTherm® Catheter is an accessory to CoreTherm® to house the Microwave Antenna and the Intraprostatic Temperature Probe for CoreTherm® treatment.

The CoreTherm® with accessories is intended for treatment of lower urinary tract symptoms (LUTS) due to benign prostatic enlargement (BPE) using microwave thermotherapy by, or under the supervision of, a trained medical practitioner.

Intended patient population

Adult male patients with BPE.

Intended user profile

To be used by, or under supervision of, a trained medical practitioner.

Intended use environment

Hospitals and health care centers.

INDICATIONS

The CoreTherm® Catheter is used for housing the Microwave Antenna and the Intraprostatic Temperature Probe to perform a CoreTherm® treatment.

See CoreTherm® IFU L-2013-009 / CoreTherm® Eagle IFU L-011, for indications to CoreTherm® treatment.

CONTRAINDICATIONS

No known contraindications with the CoreTherm® Catheter.

See CoreTherm® IFU L-2013-009 / CoreTherm® Eagle IFU L-011, for contraindications for CoreTherm® treatment.

WARNINGS

There is a risk for complications if the CoreTherm® Catheter is not handled carefully and correctly. The risk for complications can be minimized by using the product as described in this instruction for use.

Improper insertion of the CoreTherm® Catheter

Exercise care when inserting the CoreTherm® Catheter into the patient's urethra, otherwise perforation and subsequent infection of the urethra may result. If the CoreTherm® Catheter is incorrectly placed or moves during treatment, or if the Microwave Antenna comes

loose, the patient's external sphincter may overheat, causing temporary or chronic incontinence.

Use of the CoreTherm® Catheter

The CoreTherm® Catheter is intended for single use only, otherwise cross-patient urinary infection could result. Each CoreTherm® Catheter has a unique catheter serial number (SN), and the software includes a lock which prevents reuse of a catheter. Do not attempt to use a CoreTherm® Catheter for more than one treatment session.

Balloon rupture

Please make sure to check the balloon with 20 ml of air before inserting it into the patient. While inside the patient's bladder, inject 20 ml of sterile water slowly in order to prevent rupture of the balloon. Balloon rupture could lead to sliding of the catheter downwards and cause damage on the external sphincter and the urethra.

Use of saline solution

Do not use saline solution when filling the balloon and water reservoir as this may result in increased temperatures due to interaction between the microwaves and the ions.

PRECAUTIONS

- The catheter is provided sterile. Observe the expiration date on the packaging label.
- Inspect the packaging visually for any damages that could jeopardize the sterile barrier. Do not use the product if the sterile barrier has been tampered with. Damaged products shall be disposed of.

DIRECTIONS FOR USE

Preparation


CoreTherm® Catheter is supplied sterile in a sterile packaging and is a single use device.

1. Observe the serial number label on the white outer box. Save this serial number as it must later be entered before starting treatment.
2. Have your assistant offer the sealed, sterile package containing the CoreTherm® Catheter.

Check that the package is undamaged. If its sterility is compromised in any way, discard the package


and continue the session with a new, sterile CoreTherm® Catheter.

3. Handle the CoreTherm® Catheter using aseptic technique. Remove the catheter from its packaging.
4. Examine the CoreTherm® Catheter. In particular, check the balloon for leakage by inflating it with 20 ml of air. If there is any sign of damage, discard the CoreTherm® Catheter and continue the session with a new, sterile catheter. When the catheter tip is directed upwards, the outlet for the Intraprostatic Temperature Probe shall be in the 1 o'clock position.
5. Insert the disinfected Microwave Antenna into the CoreTherm® Catheter and secure the antenna by turning the Luer-lock in the clockwise direction.

NOTE: Do not lubricate the antenna.

6. Insert the sterile Intraprostatic Temperature Probe fully into the CoreTherm® Catheter, and then withdraw the probe so that its tip is visible but not protruding from the catheter.

NOTE: Do not secure the Intraprostatic Temperature Probe to the CoreTherm® Catheter yet.

NOTE: The CoreTherm® Catheter, the Microwave Antenna and the Intraprostatic Temperature Probe must be of the same version (Short or Standard) to be able to start a treatment.

7. Fill the water reservoir with 75 ml of sterile water.

Insertion of the CoreTherm® Catheter

1. Administer local anaesthetic gel into the urethra

2. Introduce the CoreTherm® Catheter into the urethra, with the tip of the catheter directed upwards all the way into the bladder in a similar manner to the insertion of a standard catheter.

NOTE: It is important that the tip of the catheter is directed upwards when the catheter is inserted into the bladder. This will ensure that the Intraprostatic Temperature Probe is inserted correctly into the left lobe of the prostate.

3. When the CoreTherm® Catheter is fully inserted, fill the balloon slowly with 20 ml of sterile water from a sterile syringe. Then it is important to check if the balloon is freely movable inside the bladder.


If the patient had a former TURP or TUMT, a special work up must be done prior to treatment. Urodynamic examination to verify obstruction and a cystoscopy to examine the bladder neck to exclude an open cavity in the bladder neck must be done. If a treatment is decided the CoreTherm® Catheter balloon is filled with 30 – 40 ml of sterile water and its correct position is always verified by TRUS before treatment start.

⚠ WARNING! Do not use saline solution when filling the balloon and water reservoir as this may result in increased temperatures due to interaction between the microwaves and the ions. Use only sterile water.

4. Gently withdraw the CoreTherm® Catheter until the balloon meets the bladder neck and is thus secured in place.

5. Hold the CoreTherm® Catheter to straighten it and insert the Intraprostatic Temperature Probe into the prostate. Advise the patient that he may experience a slight pricking sensation.
6. When fully inserted, secure the Intraprostatic Temperature Probe to the CoreTherm® Catheter by turning the locking nut of the probe connector in the clockwise direction.

An ultrasound examination may be performed to verify the position of the CoreTherm® Catheter (balloon) and the Intraprostatic Temperature Probe.

NOTE: In order to ensure correct temperature reading, check the position of the catheter regularly during the treatment by pulling it softly, in order to make sure that the catheter has not migrated.

- Place the water reservoir of the CoreTherm® Catheter into the heat exchange plate on the CoreTherm® device.

NOTE: Ensure that the water reservoir is placed flat on the heat exchange plate.

See further instructions in CoreTherm® IFU L-2013-009 / CoreTherm® Eagle IFU L-011.

Treatment

Follow CoreTherm® IFU L-2013-009 / CoreTherm® Eagle IFU L-011, on how to perform treatment.

After treatment

- Release the Intraprostatic Temperature Probe from the CoreTherm® Catheter by turning its locking nut in the counterclockwise direction.
- Carefully withdraw the Intraprostatic Temperature Probe until it entirely resides in the CoreTherm® Catheter (minimum 50 mm).

⚠ WARNING! Insufficient retraction of the Intraprostatic Temperature Probe before the catheter is removed may cause patient injuries. It is therefore important to retract the Intraprostatic Temperature Probe at least 50 mm before the catheter is removed.

- Deflate the balloon of the CoreTherm® Catheter by aspirating the fluid with a Luer syringe.
- Gently remove the catheter from the patient's urethra. When the catheter is withdrawn, it is quite normal for the patient to experience some bleeding from the urethra, due to abrasions caused by the

Intraprostatic Temperature Probe. No particular measures are necessary. The bleeding should cease spontaneously after approximately 10 minutes.

- Unlock the Microwave Antenna from the CoreTherm® Catheter as follows:
 - Hold the Luer-lock on the antenna in one hand
 - Hold the Luer part of the CoreTherm® Catheter in the other hand
 - Gently turn the Luer-lock on the antenna in the counterclockwise direction
 - When the Luer-lock is open, withdraw the Microwave Antenna from the CoreTherm® Catheter
- Make sure the Microwave Antenna is completely removed from the CoreTherm® Catheter.
- CoreTherm® Catheter is for single use only and should be disposed after use. Make sure the Microwave Antenna and the Intraprostatic Temperature Probe are completely removed from the CoreTherm® Catheter before disposal.

STORAGE

Store the CoreTherm® Catheter in a clean and dry place, with temperature between 10°C to 30°C at 10 – 80% atmospheric humidity. Expire date is labeled on the package.

SINGLE USE

CoreTherm® Catheter is for single use only.

SAFE DISPOSAL

CoreTherm® Catheter shall be disposed according to the local regulations applying to hazardous biological substances.

REPORTABILITY

Any serious incident that has occurred in relation to the medical device shall be reported to ProstaLund AB and the authority having jurisdiction in your locale.

Shaft size	22 Fr (7.3 mm)
Balloon fill	20 ml sterile water
AN connector	Luer
IP connector	Luer
Balloon inflation valve	Luer

SYMBOLS

	Manufacturer		Range of temperature
	Read instructions for use		Medical Device - Indicates the item is a medical device
	Sterilized by Ethylene Oxid		Type BF applied part
	Article number		Single sterile barrier system with protective packaging inside
	Batch number		Indicates the date when the medical device was manufactured, YYYY-MM
	Serial number		Date of expiry, YYYY-MM-DD
	Do not use if sterile barrier system or its packaging is compromised		Do not reuse. Single use device
	Range of humidity		

ORDERING INFORMATION

Device	Article Number
CoreTherm Catheter, Short	CA804120
CoreTherm Catheter, Standard	CA804220
Microwave Antenna, Short	AN806101
Microwave Antenna, Standard	AN806001
Intraprostatic Temperature Probe, AC Short (IP Short autoclavable)	IP807115
Intraprostatic Temperature Probe, AC Standard (IP Standard autoclavable)	IP807015

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